







#### The problem

Training on a manikin that has an athletes body will not prepare students for real life. Performing medical procedures could be more difficult on a robust patient.




Hyper-realistic silicone overlays for medical teams in hospitals to train. With our overlays, a greater emotional involvement is achieved, thus improving the experience of medical simulation.





Fits over your simulator with a strap to secure it in place so the simulator can be seated/ moved.





#### **Key Features**

- Corpulent body challenge
- · Works on most simulators.
- Quick and easy setup and removal.
- Medical grade silicone.
- Easy to store and maintain.
- Made with technology that prevents tears.
- Hand painted by our artists.
- Hand applied hair.
- Includes ECG connectors
- Allows auscultation
- Fits over your simulator with a strap to secure it in place so the simulator can be seated/ moved.
- For use with a simulated Defibrillator.
- Included a sustainable box for storage with support for the overlay.





After



Before



After



#### Product Codes and descriptions

The product code for the overlays is their given name followed by the selected skin tone (Light, Medium, Dark). Any additional customization will be billed separately.

Compatible with SimMan 3G, SimMan
Essential y SimMan ALS:
SimMan 3G Plus, Resusci Anne Simulator.
Juno, Ares, Apolo\*
Leonardo
\*Ask our team about limitations on Apolo



Bob



**Esther** Senior female



middle age female



**Paul** Middle age male



**Dorothy**Senior female

Compatible with Gaumard,s HAL S5301, HAL S3201 and Hal S1000.



**Esther G**Senior female

# Coming soon!

**Eva G** middle age female

# Coming soon!

**Paul G**Middle age male

# Coming soon!

**Bob G** Senior Male









How to install our overlays video



Sustainable box for storage video